



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**neolab® slide magazine with 100-200 places,  
white**

**€15.30  
plus VAT &  
Shipping**

## Product Images



## Description

---

Stackable magazine with transparent hinged lid for upright storage of slides with 76 x 26 mm. The removable insert has two rows with 50 slots each (numbered). The slots are so wide that two slides can also be inserted one behind the other (=200 slides). Without the insert, up to 400 slides can be archived in a space-saving manner. The upright storage of the slides enables them to be found quickly.

Scope of delivery: white magazine and insert, with transparent hinged lid and two index cards for filling in. Dimensions: 245 x 82 x 86 mm (L x W x H)

## Additional Information

---

No.	7-4573
Manufacturer (Brand)	neoLabLine
EAN	4058072040703
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	White
Material	Polystyrene (PS)
Wide	82 mm
Height	86 mm
Length	245 mm
Depth	245 mm
Type closure/ lid	Hinged latch/lid
suitable for	Microscope slides (76 x 26mm)
Boxes properties	stackable/collapsible with expandable capacity with labeling field coded
Type Coding	numbered
for slide dimensions	76 x 26 mm

