




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® slide magazine with 50 places, white

€11.90
plus VAT &
Shipping

Product Images

Description

Stackable magazine with 50 slots for standard slides, consisting of a box with transparent lid and removable carrier with 50 numbered slots for upright storage.

2 slides also fit in each slot.

Storage of up to 200 slides possible without carrier.

Not autoclavable.

Dimensions: 140 x 82 x 86 mm (L x W x H)

Additional Information

No.	7-7583
Manufacturer (Brand)	neoLabLine
EAN	4058072041649
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	White
Material	Polystyrene (PS)
Wide	82 mm
Height	86 mm
Length	140 mm
Depth	140 mm
Type closure/ lid	Hinged latch/lid
suitable for	Microscope slides (76 x 26mm)
Boxes properties	stackable/pluggable with labeling field coded
Type Coding	numbered
for slide dimensions	76 x 26 mm

