



neoLab® Slide Folders, white

Product Images



Description

Slide folders for safe and clean storage or transport of 76 x 26 mm slides.
With 20 numbered spaces.
Transparent, hinged front for easy identification of slides from the outside.
Easy removal of slides, stackable.

No.	7-4597
Manufacturer (Brand)	neoLab
EAN	4058072040895
Color	Weiß
Material	Polystyrol (PS)
Wide	192 mm
Height	11 mm
Length	295 mm
Depth	295 mm
Type closure/ lid	Push-on/Steckverschluss
suitable for	Objektträger (76 x 26mm)
Boxes properties	stapel-/zusammensteckbar für den Versand geeignet codiert
Type Coding	nummeriert
for slide dimensions	76 x 26 mm



Additional Information