



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neolab® slide shipping container green, 50
pcs./pack**

**€30.40
plus VAT &
Shipping**

Product Images



Description

Green slide holder for 1 slide, can be used as a shipping container or also as a secure stacking base.

Easy and safe to close, easy removal of the slide, no scratching of the surface.

A special mechanism lets stacked holders snap into each other.

Suitable for slides 76 x 26 mm, also with cover slip.

Temperature resistant from -80°C to +80°C.

Outer dimensions: 89 x 29 x 6 mm (L x W x H)

PU: 50 pcs. /pack

Additional Information

No.	2-1825
Manufacturer (Brand)	neoLab
EAN	4058072058074
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	green
Material	Polystyrene (PS)
Wide	29 mm
Height	6 mm
Length	89 mm
Depth	89 mm
Type closure/ lid	Push-on/plug-in closure
suitable for	Microscope slides (76 x 26mm)
Boxes properties	stackable/collapsible with expandable capacity suitable for shipping
for slide dimensions	76 x 26 mm

