



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Tipping box set made of PC, 6 containers

**Price on
request
plus VAT &
Shipping**

Product Images



Description

The neoLab tilting boxes provide a full overview at all times

The contents can be seen immediately due to the transparent box made of robust polycarbonate. Complete flexibility is offered when filling or removing the contents: The boxes can be easily tilted forward or removed completely.

- All versions with a width of 60 cm
- Particularly scratch-resistant
- Individually stackable

Additional Information

No.	2-1651
Manufacturer (Brand)	neoLab
EAN	4058072156602
Transport temperature	Room temperature
Color	transparent (crystal clear)
Material	Polycarbonate (PC)
Wide	91 mm
Height	113 mm
Length	600 mm
Depth	600 mm
Boxes properties	stackable / interlockable
Number of drawers	6

