



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® Tipping box set made of PC, 3 containers

**€79.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

### The neoLab tilting boxes provide a full overview at all times

The contents can be seen immediately due to the transparent box made of robust polycarbonate. Complete flexibility is offered when filling or removing the contents: The boxes can be easily tilted forward or removed completely.

- All versions with a width of 60 cm
- Particularly scratch-resistant
- Individually stackable

## Additional Information

---

No.	2-1654
Manufacturer (Brand)	neoLab
EAN	4058072157203
Transport temperature	Room temperature
Color	transparent (crystal clear)
Material	Polycarbonate (PC)
Wide	197 mm
Height	240 mm
Length	600 mm
Depth	600 mm
Boxes properties	stackable / interlockable
Number of drawers	3

