



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Small parts magazine with 24 drawers 175 x 55 x 37 mm

€69.70
plus VAT &
Shipping

Product Images



Description

These magazines can be placed individually or attached to the wall.
With sturdy steel frame, powder coated light gray.
Drawers made of crystal clear polystyrene with large handle.

External dimensions of the magazines

Order No. 2-1311

16 x 30 x 27 cm (D x W x H)

Additional Information

No.	2-1311
Manufacturer (Brand)	neolab
EAN	4058072056018
Transport temperature	Room temperature
Color	gray transparent (crystal clear)
Material	Polystyrene (PS) steel powder coated
Wide	160 mm
Height	300 mm
Length	270 mm
Depth	270 mm
suitable for	variable applications
Boxes properties	stackable / interlockable
Number of drawers	30
Width drawer	175.0 mm
Height drawer	37.0 mm
Length drawer	55.0 mm

