



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® Small parts magazine with 4 drawers 175 x 277 x 59 mm

**€68.50**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

These magazines can be placed individually or attached to the wall.  
With sturdy steel frame, powder coated light gray.  
Drawers made of crystal clear polystyrene with large handle.

External dimensions of the magazines

Order No. 2-1315

16 x 30 x 27 cm (D x W x H)

## Additional Information

---

No.	2-1315
Manufacturer (Brand)	neoLab
EAN	4058072056056
Transport temperature	Room temperature
Color	gray transparent (crystal clear)
Material	Polystyrene (PS) steel powder coated
Wide	160 mm
Height	300 mm
Length	270 mm
Depth	270 mm
suitable for	variable applications
Boxes properties	stackable / interlockable
Number of drawers	4
Width drawer	175.0 mm
Height drawer	59.0 mm
Length drawer	277.0 mm

