




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neolab® small parts magazine with 60 drawers
175 x 55 x 37 mm**

**€116.10
plus VAT &
Shipping**



Description

These magazines can be placed individually or attached to the wall.
With sturdy steel frame, powder coated light gray.
Drawers made of crystal clear polystyrene with large handle.

Outside dimensions of the magazines

Order No. 2-1611

16 x 30 x 54 cm (D x W x H)

Additional Information

| | |
|-----------------------|--------------------------------------|
| No. | 2-1611 |
| Manufacturer (Brand) | neoLab |
| EAN | 4058072057664 |
| Transport temperature | Room temperature |
| Color | gray transparent (crystal clear) |
| Material | Polystyrene (PS) steel powder coated |
| Wide | 160 mm |
| Height | 300 mm |
| Length | 540 mm |
| Depth | 540 mm |
| suitable for | variable applications |
| Boxes properties | stackable / interlockable |
| Number of drawers | 60 |
| Width drawer | 175.0 mm |
| Height drawer | 37.0 mm |
| Length drawer | 55.0 mm |

