



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® small parts magazine with 24 drawers 175 x 92 x 59 mm

€115.80
plus VAT &
Shipping

Product Images



Description

These magazines can be placed individually or attached to the wall.
With sturdy steel frame, powder coated light gray.
Drawers made of crystal clear polystyrene with large handle.

Outside dimensions of the magazines

Order No. 2-1614

16 x 30 x 54 cm (D x W x H)

Additional Information

| | |
|-----------------------|--------------------------------------|
| No. | 2-1614 |
| Manufacturer (Brand) | neoLab |
| EAN | 4058072057695 |
| Transport temperature | Room temperature |
| Color | gray transparent (crystal clear) |
| Material | Polystyrene (PS) steel powder coated |
| Wide | 160 mm |
| Height | 300 mm |
| Length | 540 mm |
| Depth | 540 mm |
| suitable for | variable applications |
| Boxes properties | stackable / interlockable |
| Number of drawers | 24 |
| Width drawer | 175.0 mm |
| Height drawer | 59.0 mm |
| Length drawer | 92.0 mm |

