



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Magazine for embedding cassettes, PS

€59.00
plus VAT &
Shipping

Product Images



Description

Stackable drawers for storing 250 embedding cassettes.

With labeling labels.

Material: PS

Dimensions: 405 x 230 x 51 mm (L x W x H)

Additional Information

No.	7-0020
Manufacturer (Brand)	neoLab
EAN	4058072026592
Transport temperature	Room temperature
Color	Black
Material	Polystyrene (PS)
Wide	230 mm
Height	51 mm
Length	405 mm
Depth	405 mm
suitable for	Embedding cassettes
Capacity	250

