



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Assortment cabinet with 6 drawers w. compartment division

€148.70
plus VAT &
Shipping

Product Images



Description

Order No. 2-1137 Assortment cabinet

6drawers:

2 x 30 compartments with 44 x 50 mm

2 x 16 compartments with 55 x 77 mm

2 x 12 compartments with 71 x 75 mm

Internal compartment height: 40 mm

Case: anthracite, stackable

Basic dimensions: 265 x 305 x 355 mm

Drawers: assorted colors

Material: PS

Additional Information

| | |
|-----------------------|---|
| No. | 2-1137 |
| Manufacturer (Brand) | neoLab |
| EAN | 4058072055608 |
| Transport temperature | Room temperature |
| Color | anthracite blue green red yellow grey white |
| Material | Polystyrene (PS) |
| Wide | 355 mm |
| Height | 265 mm |
| Length | 305 mm |
| Depth | 305 mm |
| suitable for | variable applications |
| Boxes properties | stackable/assemblable with removable insert with insert with labeling field |
| Number of drawers | 6 |
| Width drawer | 350.0 mm |
| Height drawer | 40.0 mm |
| Length drawer | 300.0 mm |
| Type insert | with compartments |

