



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Assortment cabinet with 4 drawers with compartment division, high

€170.30
plus VAT &
Shipping

Product Images



Description

Order No. 2-1138 Assortment cabinet

4drawers:

2x 6 compartments with 75 x 160 mm

2x 12 compartments with 71 x 75 mm

Internal compartment height: 62 mm

Case: anthracite, stackable

Basic dimensions: 265 x 305 x 355 mm

Drawers: assorted colors

Material: PS

Additional Information

No.	2-1138
Manufacturer (Brand)	neoLab
EAN	4058072055615
Transport temperature	Room temperature
Color	anthracite blue green red yellow grey white
Material	Polystyrene (PS)
Wide	355 mm
Height	265 mm
Length	305 mm
Depth	305 mm
suitable for	variable applications
Boxes properties	stackable/assemblable with removable insert with insert with labeling field
Number of drawers	4
Width drawer	350.0 mm
Height drawer	62.0 mm
Length drawer	300.0 mm
Type insert	with compartments

