



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Alu basket medium, 256 x 152 x 152 mm

€44.70
plus VAT &
Shipping

Product Images



Description

Ideal for rinsing and storing test tubes!

Lightweight but particularly durable aluminum baskets for many laboratory applications.
Can be used in dishwashers and autoclaves.

Additional Information

No.	2-5091
Manufacturer (Brand)	neoLab
EAN	4058072004309
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	silver
Material	Aluminum
Volume	5.9 l
Wide	152 mm
Form	rectangular
Height	152 mm
Length	256 mm
Depth	256 mm
TBST MAX	121 °C
TBST MIN	-20 °C
Size (small, medium, large)	medium

