



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® plastic basket 127 x 102 x 102 mm, with lid

**€46.50**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

For storing, rinsing, autoclaving or incubating e.g. reagent or culture tubes, Petri dishes etc.

Material: polypropylene, autoclavable up to 121°C

Very stable, high chemical resistance, can be stored stacked inside each other.

## Additional Information

---

No.	1-3123
Manufacturer (Brand)	neoLab
EAN	4058072035143
autoclavable at 121°C	Yes
Color	White
Material	Polypropylene (PP)
Volume	1.3 l
Wide	127 mm
Form	rectangular
Height	102 mm
Length	102 mm
Depth	102 mm
Baskets features	with special chemical resistance with closure/lid
TBST MAX	121 °C
TBST MIN	0 °C

