



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® plastic basket 152 x 152 x 152 mm, with lid

€65.00
plus VAT &
Shipping

Product Images



Description

For storing, rinsing, autoclaving or incubating e.g. reagent or culture tubes, Petri dishes etc.

Material: polypropylene, autoclavable up to 121°C

Very stable, high chemical resistance, can be stacked inside each other for storage.

Additional Information

No.	1-3124
Manufacturer (Brand)	neoLab
EAN	4058072035150
autoclavable at 121°C	Yes
Color	White
Material	Polypropylene (PP)
Volume	3.5 l
Wide	152 mm
Form	square
Height	152 mm
Length	152 mm
Depth	152 mm
TBST MAX	121 °C
TBST MIN	0 °C
Baskets features	with special chemical resistance with closure/lid

