



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® plastic basket 229 x 229 x 229 mm, with lid

**€53.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

For storing, rinsing, autoclaving or incubating e.g. reagent or culture tubes, Petri dishes etc.

Material: polypropylene, autoclavable up to 121°C

Very stable, high chemical resistance, can be stored stacked inside each other.

## Additional Information

---

|                       |   |
|-----------------------|---|
| No.                   | 1-3125  |
| Manufacturer (Brand)  | neoLab  |
| EAN                   | 4058072035167                                     |
| autoclavable at 121°C | Yes   |
| Color                 | White   |
| Material              | Polypropylene (PP)                                |
| Volume                | 12 l  |
| Wide                  | 229 mm  |
| Form                  | square  |
| Height                | 229 mm  |
| Length                | 229 mm  |
| Depth                 | 229 mm  |
| TBST MAX              | 121 °C  |
| TBST MIN              | 0 °C  |
| Baskets features      | with special chemical resistance with closure/lid |

