



neoLab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e


Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Left stainless steel wire basket, 12 x 12 x 12 cm

**€36.00** As low  
as  
**€34.00**  
plus VAT &  
Shipping

### Product Images



## Description

---

Fine mesh stainless steel baskets for storage and transport. The 5 mm thick frame provides high stability, the 1 mm thick wire mesh can withstand quite a bit!

Material: 4301 stainless steel, finely electropolished

Mesh size: 8 mm

## Additional Information

---

No.	LR-0013
Manufacturer (Brand)	Linker
old neolab article no.	E-7020
EAN	4058072288310
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	silver
Material	Stainless steel
Volume	1.7 l
Wide	120 mm
Form	square
Height	120 mm
Length	120 mm
Depth	120 mm
TBST MAX	121 °C
TBST MIN	-20 °C
Baskets features	with special chemical resistance

