



neoLab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Left stainless steel wire basket, 14 x 14 x 14 cm

€41.00 As low
as
€37.00
plus VAT &
Shipping

Product Images



Description

Fine mesh stainless steel baskets for storage and transport. The 5 mm thick frame provides high stability, the 1 mm thick wire mesh can withstand quite a bit!

Material: 4301 stainless steel, finely electropolished

Mesh size: 8 mm

Additional Information

No.	LR-0014
Manufacturer (Brand)	Linker
old neolab article no.	E-7021
EAN	4058072288327
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	silver
Material	Stainless steel
Volume	2.7 l
Wide	140 mm
Form	square
Height	140 mm
Length	140 mm
Depth	140 mm
TBST MAX	121 °C
TBST MIN	-20 °C
Baskets features	with special chemical resistance

