



neoLab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Left stainless steel wire basket, 16 x 16 x 16 cm

€47.00 As low

as

€44.00

**plus VAT &
Shipping**

Product Images



Description

Fine mesh stainless steel baskets for storage and transport. The 5 mm thick frame provides high stability, the 1 mm thick wire mesh can withstand quite a bit!

Material: 4301 stainless steel, finely electropolished

Mesh size: 8 mm

Additional Information

No.	LR-0015
Manufacturer (Brand)	Linker
old neolab article no.	E-7022
EAN	4058072288334
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	silver
Material	Stainless steel
Volume	3.8 l
Wide	160 mm
Form	square
Height	160 mm
Length	160 mm
Depth	160 mm
TBST MAX	121 °C
TBST MIN	-20 °C
Baskets features	with special chemical resistance

