



neoLab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



RSG stainless steel wire basket, 30 x 20 x 20 cm

€90.57
plus VAT &
Shipping

Product Images



Description

Fine mesh stainless steel baskets for storage and transport. The 5 mm thick frame provides high stability, the 1 mm thick wire mesh can withstand quite a bit!

Material: 4301 stainless steel, finely electropolished

Mesh size: 8 mm

Additional Information

No.	RG-1029
Manufacturer (Brand)	RSG
old neolab article no.	E-7025
EAN	4058072338619
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	silver
Material	Stainless steel
Volume	12 l
Wide	200 mm
Form	rectangular
Height	200 mm
Length	300 mm
Depth	300 mm
TBST MAX	121 °C
TBST MIN	-20 °C
Baskets features	with special chemical resistance

