



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



DURAN® Bottle carrier system GL 45, PP, 5.0 L, yellow

€119.73
plus VAT &
Shipping

Product Images



Description

The system is based on a GL 45 bottle adapter specially designed to securely fasten the sling system around the neck of the bottle. Strong, resistant straps and handles made of non-slip material always hold the bottle securely in the desired position. Further protection is provided by the shock-absorbing silicone carrying base. All components of the bottle carrying system are autoclavable and can be easily cleaned, even when attached to the bottle.

The DURAN® GL 45 bottle carrying system is available for the following four bottle sizes: 2, 5, 10 and 20 liters and comes in three different designs depending on the bottle size used.

The **5-liter carrying system**, for example, allows safe and controlled pouring thanks to an additional pouring handle on the side.

Additional Information

No.	SD-2804
Manufacturer (Brand)	DURAN
EAN	4058072364038
Transport temperature	Room temperature
Color	Black Yellow
Material	Polypropylene (PP)
Volume	5 l
Form	Strips

