



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® Bottle Carrying Rack (PP) for 6 bottles up to 500 ml

**€388.40**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Carrying rack made of polypropylene for bottles.

The basic rack is without division and can therefore also be used for other transport tasks.

Round, comfortable carrying handle with 25 mm Ø facilitates transport.

Inner dimensions: 36 x 18 x 15 cm (L x W x H)

Accessories: Insert with 6 compartments (Order No. 2-1128)

## Additional Information

---

No.	2-1127
Manufacturer (Brand)	neoLab
EAN	4058072055578
Color	White
Material	Low-density polyethylene (LDPE)
Wide	180 mm
Width/inside	155 mm
Form	angular
Height	150 mm
Length	360 mm
Depth	360 mm
suitable for	Bottles. square Bottles. round
Length inside	335 mm
for vessels with volume up to	500 ml
for vessels with DM up to	80.0 mm
for type bottle	Laboratory bottles. round square bottles

