



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**neolab® Safety box made of EPS for 1 x 100 ml,  
100 pcs./pack**

**€147.90  
plus VAT &  
Shipping**

## Product Images



## Description

---

Safety boxes made of polystyrene (EPS) for transporting square and round glass bottles.

Flame retardant.

Optimum utilization of the EPS cellular framework, thus high strength.

No absorption of moisture.

Light weight saves transport costs.

Heat resistance up to +80°C.

**Order No. 2-2801**

**Capacity:** 1 bottle

## Additional Information

---

No.	2-2801
Manufacturer (Brand)	neoLab
EAN	4058072069155
Transport temperature	Room temperature
Color	White
Material	Polystyrene (EPS. Styropor®)
Wide	80 mm
Form	angular
Height	155 mm
Length	80 mm
Depth	80 mm
suitable for	Bottles. round Bottles. square
Baskets features	with closure/lid
for vessels with volume up to	100 ml
for vessels with DM up to	50.0 mm

