



neoLab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**neoLab® Safety box made of EPS for 1 x 250 ml,  
54 pcs./pack**

**€112.90  
plus VAT &  
Shipping**

## Product Images



## Description

---

Safety boxes made of polystyrene (EPS) for transporting square and round glass bottles.

Flame retardant.

Optimum utilization of the EPS cellular framework, thus high strength.

No absorption of moisture.

Light weight saves transport costs.

Heat resistance up to +80°C.

**Order No. 2-2802**

**Capacity:** 1 bottle

## Additional Information

---

No.	2-2802
Manufacturer (Brand)	neoLab
EAN	4058072069162
Transport temperature	Room temperature
Color	White
Material	Polystyrene (EPS. Styropor®)
Wide	100 mm
Form	angular
Height	190 mm
Length	100 mm
Depth	100 mm
suitable for	Bottles. round Bottles. square
Baskets features	with closure/lid
for vessels with volume up to	250 ml
for vessels with DM up to	70.0 mm

