



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neolab® Safety box made of EPS for 1 x 500 ml,
50 pcs./pack**

**€135.00
plus VAT &
Shipping**

Product Images



Description

Safety boxes made of polystyrene (EPS) for transporting square and round glass bottles.

Flame retardant.

Optimum utilization of the EPS cellular framework, thus high strength.

No absorption of moisture.

Light weight saves transport costs.

Heat resistance up to +80°C.

Order No. 2-2803

Capacity: 1 bottle

Additional Information

No.	2-2803
Manufacturer (Brand)	neoLab
EAN	4058072069179
Transport temperature	Room temperature
Color	White
Material	Polystyrene (EPS. Styropor®)
Wide	116 mm
Form	angular
Height	230 mm
Length	116 mm
Depth	116 mm
suitable for	Bottles. round Bottles. square
Baskets features	with closure/lid
for vessels with DM up to	84.0 mm
for vessels with volume up to	500 ml

