



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Safety box made of EPS for 1 x 2500 ml (139 mm Ø, 300 mm high)

€9.72
plus VAT &
Shipping

Product Images



Description

Safety boxes made of polystyrene (EPS) for transporting square and round glass bottles.

Flame retardant.

Optimum utilization of the EPS cellular framework, thus high strength.

No absorption of moisture.

Light weight saves transport costs.

Heat resistance up to +80°C.

Order No. 2-2805

Capacity: 1 bottle

Additional Information

No.	2-2805
Manufacturer (Brand)	neoLab
EAN	4058072069193
Transport temperature	Room temperature
Color	White
Material	Polystyrene (EPS. Styropor®)
Wide	187 mm
Form	angular
Height	380 mm
Length	187 mm
Depth	187 mm
suitable for	Bottles. round Bottles. square
Baskets features	with closure/lid
for vessels with volume up to	2500 ml
for vessels with DM up to	139.0 mm

