



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Safety box made of EPS for 6 x 100 ml

€4.35
plus VAT &
Shipping

Product Images



Description

Safety boxes made of polystyrene (EPS) for transporting square and round glass bottles.

Flame retardant.

Optimum utilization of the EPS cellular framework, thus high strength.

No absorption of moisture.

Light weight saves transport costs.

Heat resistance up to +80°C.

Order No. 2-2808

Capacity: 6 bottles

Additional Information

| | |
|-------------------------------|--------------------------------|
| No. | 2-2808 |
| Manufacturer (Brand) | neoLab |
| EAN | 4058072069209 |
| Transport temperature | Room temperature |
| Color | White |
| Material | Polystyrene (EPS. Styropor®) |
| Wide | 159 mm |
| Form | angular |
| Height | 205 mm |
| Length | 143 mm |
| Depth | 143 mm |
| suitable for | Bottles. round Bottles. square |
| Baskets features | with closure/lid |
| for vessels with volume up to | 100 ml |
| for vessels with DM up to | 50.0 mm |

