



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



DURAN® filtering beaker with graduation and spout, thick-walled, 500 ml

€338.10
plus VAT &
Shipping

Product Images



Description

DURAN® beakers, thick-walled, 500 ml Beakers are cylindrical vessels with a straight side wall. Thanks to this shape, beakers are universally applicable laboratory vessels and are usually used as boiling vessels, as collecting vessels and for preparing, mixing and heating solutions. The large opening of the beaker makes it easy to fill. Beakers are available in different designs and shapes: high and low form, with and without spout, with and without scale. Special features and advantages: Thick-walled (filtration beaker) Precise spout for clean working Due to reduced thermal shock resistance, can only be used as a boiling glass to a limited extent Easily readable scale with line graduation for measuring the contents Burnt-in and therefore permanently durable white printed image Application: Filtering, heating, diluting solutions, mixing several components, measuring liquid quantities, sizes 10000, 15000 and 20000 ml without graduation.

Additional Information

No.	SD-0244
Reg. trademark	DURAN ®
Manufacturer (Brand)	DURAN
VGKL number	112113144
EAN	4058072182830
Manufacturer EAN	4032051003477
Gross weight	0.315
Transport temperature	Room temperature
Storage temperature	Room temperature
Net weight	0.30100 kg
Material	Borosilicate glass 3.3
Volume	500 ml
sterile	No
Country of origin	Croatia
Packaging volume	0.01286300 cc
Customs tariff number	70172000

