



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Blotting Container, PP, 16 x 110 mm (with accessories) , 5 pcs/pack

€56.90
plus VAT &
Shipping

Product Images



Description

The neoLab blotting containers made of PP are the solution for optimal results from Western Blot, Northern Blot and other blotting systems.

Depending on how many lanes of blotting paper are to be placed, a container with the optimal matching dimensions can be selected, which are tuned to ensure a reduction of the required antibodies.

The containers have a lid and small knobs on the outer bottom for stacking.

Order #: 1-1344 to 1-1346

Minimize the use of expensive antibodies during incubation with an optimally sized tray.

Polypropylene trays are available for common midi and mini gels.

Stackable, resealable, reusable and autoclavable.

All trays have an attached lid to protect against spills.

Single lane trays are supplied as a group of five connected lanes; these can be easily separated into individual lanes prior to use.

Additional Information

No.	1-1344
Manufacturer (Brand)	neoLab
EAN	4058072177553
autoclavable at 121°C	Yes
Color	transparent (blue)
Material	Polypropylene (PP)
sterile	No
Wide	16 mm
Form	rectangular
Length	115 mm
Closure / lid properties	closed (attached)
Type bowl, dish, crucible	Blotting tank

