



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e


Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® Blotting Container for Midi Blots, PVC, Set of 3

**€58.20**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

The neoLab blotting containers made of transparent PVC, are the solution for optimal results from Western Blot, Northern Blot and other blotting systems.

Depending on how many lanes of blotting paper are to be placed, a container with the optimal matching dimensions can be selected, which are tuned to ensure a reduction of the required antibodies.

The containers have a lid and small knobs on the outer bottom for stacking

### **Order No.: 1-1348**

Sets of 2 sizes can be used for maximum efficiency in Western blots.

Reduce antibody volume by using only the tray you need.

Sets include full and half blot trays as well as single lane, 2 lane or 3 lane (midi only)

The included lid protects samples and allows stacking.

All 15.4 cm long

- 18.4 cm (3 pieces, incl. 3 lids)
- 8.4 cm (6 pieces)
- 5.8 cm (9 pieces)
- 3.7 cm (12 pieces)
- 2.6 cm (18 pieces)

## Additional Information

---

No.	1-1348
Manufacturer (Brand)	neoLab
EAN	4058072171186
autoclavable at 121°C	No
Color	transparent (milky)
Material	Polyvinyl chloride (PVC. non-specific)
sterile	No
Form	rectangular
Closure / lid properties	overlying
Type bowl, dish, crucible	Blotting tank

