



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Nalgene large polypropylene sterilization trays Case of 6 9.5 l

**€56.00** As low  
as  
**€45.17**  
plus VAT &  
Shipping

### Product Images

## Description

---

Nalgene™ Large Polypropylene Sterilization Trays Dry and sterilize instruments and laboratory items with the Thermo Scientific™ Nalgene™ Large Polypropylene Sterilization Trays. Description: Large sterilization trays are useful for drying and sterilizing instruments and laboratory equipment or for collecting and washing soiled laboratory items Dimensions (L x W x H): 361 x 311 x 143 mm (14.25 x 12.25 x 5.675 in.) Height (Anglo): 2.5 in. Length (Anglo): In 10.25. Quantity Per Carton: 6 Quantity Per Package: 1 Shape: Rectangular Width (Anglo): 6.25 in. Autoclavable: Autoclavable Type: Sterilization Tray Unit Size: Case of 6

## Additional Information

---

No.	CD-1791
Manufacturer (Brand)	Nalgene
VGKL number	238769001
EAN	4066292480193
Gross weight	3.72
Transport temperature	Room temperature
Material	Polypropylene (PP)
Volume	9500 ml
Country of origin	USA
Packaging volume	0.03900000 cc
Customs tariff number	39269097

