



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



DURAN® Crystallizing dish, with spout, 3500 ml

€129.70 As low
as
€99.77
plus VAT &
Shipping

Product Images



Description

DURAN® Crystallizing dish, with spout, 3500 ml Converting substances and liquids into other states of aggregation and thus initiating further reactions is very important for daily work in the laboratory. A wide variety of vessels, apparatus and equipment are used for these activities. Special features and advantages: Flat and wide shape enable fast and uniform evaporation of solvents Large base area Precise spout for clean work Burnt-in and therefore permanently durable white print Application: crystallization , preparation of concentrated solutions, Available with spout DIN 12338 and without spout DIN 12337.

Additional Information

No.	SD-0459
Reg. trademark	DURAN ®
Manufacturer (Brand)	DURAN
VGKL number	112131163
EAN	4058072175405
Manufacturer EAN	4032051008366
Gross weight	0.708
Transport temperature	Room temperature
Storage temperature	Room temperature
Net weight	0.62000 kg
Material	Borosilicate glass 3.3
Volume	3500 ml
sterile	No
Country of origin	Croatia
Packaging volume	0.06303800 cc
Customs tariff number	70172000

