



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Moonlab® Petri dish, Ø: 90 mm, height: 15 mm,
3 cams, clear, sterile, PS, individually packed,
PU:10 pcs/bag**

€5.37 As low as
€4.64
**plus VAT &
Shipping**

Product Images



Description

Disposable Petri dishes, sterile made of sturdy polystyrene. Ideally suited for all work in microbiology. High optical transparency for better examination. The lid has 3 cams for optimal air circulation. Temperature resistant up to 55 °C.

Packing unit: 10 pieces/bag (individually packed)

Additional Information

No.	4-0031
Manufacturer (Brand)	Moonlab
EAN	4058072085339
Transport temperature	Room temperature
diameter	90 mm
Material	Polystyrene (PS)
sterile	Yes
Height	15 mm
Bowls, pots and jars properties	with closure/lid with knobs RNA-free DNA-free
Type of sterility	radiation sterile
Growth area	63.62 cm ²
Format/ Packing	individually packed

