



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neoCulture cell culture dishes 35 mm Ø, PS, sterile, box of 960 pieces (96 x 10 pcs)

Special Price

€119.00 was

€215.00

**plus VAT &
Shipping**

Product Images



Description

Cell and tissue culture dishes with treated surface are ideal for adherent cell growth. They can also be used for sample separation, pretreatment, storage and so on. The optically clear surface of the cell culture dishes enables distortion-free microscopy.

- Smooth, streak-free dish surface
- Easily stackable
- Lid with cams for effective gas exchange
- Suitable for many common vending machines
- Sterilized by gamma radiation
- DNase- / RNase-free
- Pyrogen-free

Additional Information

No.	C-8211_RELAB
neoGreen	Yes
Manufacturer (Brand)	re:Lab
EAN	4066292691018
Transport temperature	Room temperature
diameter	35 mm
Material	Polystyrene (PS)
sterile	Yes
Height	12.5 mm
Bowls, pots and jars properties	surface optimized
Growth area	8.50 cm ²

