



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Lenz glass filter plate, Ø 40 mm, por. 2

**€56.10**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Glass filter plates

Filter plates, centered, rim unfused, material: borosilicate glass 3.3

All glass filters are made of pure and high quality borosilicate glass 3.3 without added binders and can be fused directly with all known laboratory glasses of this quality.

## Additional Information

---

No.	LZ-2246
Manufacturer (Brand)	Lenz
old neoLab article no.	B-4250
EAN	4058072319793
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	White
Material	Borosilicate glass
DM inside	40 mm
Floor outside	flat
Porosity	2
Type filter	Glass fiber filter

