



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Lenz glass filter plate, Ø 30 mm, por. 4

€50.50
plus VAT &
Shipping

Product Images



Description

Glass filter plates

Filter plates, centered, rim unfused, material: borosilicate glass 3.3

All glass filters are made of pure and high quality borosilicate glass 3.3 without added binders and can be fused directly with all known laboratory glasses of this quality.

Additional Information

No.	LZ-2257
Manufacturer (Brand)	Lenz
old neoLab article no.	B-4261
EAN	4058072319908
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	White
Material	Borosilicate glass
DM inside	30 mm
Floor outside	flat
Porosity	4
Type filter	Glass fiber filter

