



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neolab® stainless steel pot 4 l, with lid, 20 cm
inner Ø, 15 cm high**

**€56.80
plus VAT &
Shipping**

Product Images



Description

18/10 stainless steel, heavy quality, brushed satin finish with high-gloss rim.
Encapsulated sandwich base, suitable for induction.
With tight-fitting lid.
Solid cold handles.

Additional Information

No.	E-1123
Manufacturer (Brand)	neolab
EAN	4058072063023
Transport temperature	Room temperature
diameter	240 mm
Color	silver
Material	Stainless steel
Volume	7 l
Height	170 mm
suitable for	Induction stoves
Bowls, pots and jars properties	with closure/lid
Scope of delivery	Stainless steel pot, matching lid

