



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



DURAN® watch glass bowl, rim fused, d = 80 mm

€17.03 As low as
€13.10
plus VAT &
Shipping

Product Images



Description

DURAN® watch glass dish, Ø 80 mm The multitude of different work steps in the laboratory environment demands a variety of universally applicable glassware. These glassware devices usually have a fairly simple shape and are characterized by their versatility. Watch glasses, for example, are concavely curved. The small trough thus created allows liquids or solids to be poured in. Special features and advantages: Reduced risk of injury due to fused rim
Application: Weighing, covering, pre-sampling, drying of materials, Available in Boro 3.3 or soda-lime glass

Additional Information

No.	SD-0474
Reg. trademark	DURAN ®
Manufacturer (Brand)	DURAN
VGKL number	112132141
EAN	4058072181499
Manufacturer EAN	4032051008571
Gross weight	0.026
Storage temperature	Room temperature
Transport temperature	Room temperature
Net weight	0.02400 kg
diameter	80 mm
Material	Borosilicate glass 3.3
Country of origin	Germany
sterile	No
Region of origin	Baden-Württemberg
Customs tariff number	70172000

