



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



DURAN® watch glass bowl, rim fused, d=100mm, soda-lime glass

€4.45 As low as
€3.42
plus VAT &
Shipping

Product Images

Description

Watch glass dish, soda-lime glass, Ø 100 mm The multitude of different work steps in the laboratory environment demands a variety of universally applicable glassware. These glassware devices usually have a rather simple shape and are characterized by their versatility. Watch glasses, for example, are concavely curved. The small trough thus created allows liquids or solids to be poured in. Special features and advantages: Reduced risk of injury due to fused rim Application: Weighing, covering, pre-sampling, drying of materials, Available in Boro 3.3 or soda-lime glass

Additional Information

No.	SD-0989
Reg. trademark	DURAN ®
Manufacturer (Brand)	DURAN
VGKL number	112332146
EAN	4058072165659
Manufacturer EAN	4032051023413
Gross weight	0.04
Storage temperature	Room temperature
Transport temperature	Room temperature
Net weight	0.03800 kg
diameter	100 mm
Material	Lime soda glass
Country of origin	Czech Republic
sterile	No
Packaging volume	0.00046700 cc
Customs tariff number	70179000

