



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



DURAN® watch glass bowl, rim fused, d=40mm, soda-lime glass

€23.10
plus VAT &
Shipping

Product Images



Description

Watch glass dish, soda-lime glass, Ø 40 mm The multitude of different work steps in the laboratory environment demands a variety of universally applicable glassware. These glassware devices usually have a rather simple shape and are characterized by their versatility. Watch glasses, for example, are concavely curved. The small trough created in this way allows liquids or solids to be poured in. Special features and advantages: Reduced risk of injury due to fused rim Application: Weighing, covering, pre-sampling, drying of materials, Available in Boro 3.3 or soda-lime glass

Additional Information

No.	SD-0983
Reg. trademark	DURAN ®
Manufacturer (Brand)	DURAN
VGKL number	112332124
EAN	4058072164577
Manufacturer EAN	4032051023352
Gross weight	0.005
Transport temperature	Room temperature
Storage temperature	Room temperature
diameter	40 mm
Net weight	0.00400 kg
Material	Lime soda glass
sterile	No
Country of origin	Czech Republic
Packaging volume	0.00008100 cc
Customs tariff number	70179000

