



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## DURAN® watch glass bowl, rim fused, d=70mm, soda-lime glass

**€28.30**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Watch glass bowl, Klak soda glass, Ø 70 mm The multitude of different work steps in the laboratory environment demands a variety of universally applicable glassware. These glassware devices usually have a rather simple shape and are characterized by their versatility. Watch glasses, for example, are concavely curved. The small trough thus created allows liquids or solids to be poured in. Special features and advantages: Reduced risk of injury due to fused rim Application: Weighing, covering, pre-sampling, drying of materials, Available in Boro 3.3 or soda-lime glass

## Additional Information

---

No.	SD-0986
Reg. trademark	DURAN ®
Manufacturer (Brand)	DURAN
VGKL number	112332138
EAN	4058072182960
Manufacturer EAN	4032051023383
Gross weight	0.022
Storage temperature	Room temperature
Transport temperature	Room temperature
Net weight	0.02000 kg
diameter	70 mm
Material	Lime soda glass
Country of origin	Czech Republic
sterile	No
Packaging volume	0.00028800 cc
Customs tariff number	70179000

