



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Lenz weighing glass with NS lid stopper, Ø 30 x
H. 50 mm, ground joint NS 29/12, high form**

**€97.68
plus VAT &
Shipping**

Product Images

Description

Weighing glasses with NS lid stopper

High form

With replaceable NS lid stopper, with labeling field. Dimensions (mm) and contents (ml) permanently branded.

Additional Information

No.	LZ-1515
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3518
EAN	4058072312541
Transport temperature	Room temperature
diameter	30 mm
Color	transparent (clear)
Material	Borosilicate glass
Volume	20 ml
Form	high shape
Height	50 mm
Type closure/ lid	NS lid plug

