



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Lenz weighing glass with NS lid stopper, Ø 35 x
H. 50 mm, ground joint NS 34/12, high form**

**€104.04
plus VAT &
Shipping**



Description

Weighing glasses with NS lid stopper

High form

With replaceable NS lid stopper, with labeling field. Dimensions (mm) and contents (ml) permanently branded.

Additional Information

| | |
|------------------------|---------------------|
| No. | LZ-1516 |
| Manufacturer (Brand) | Lenz |
| old neoLab article no. | B-3519 |
| EAN | 4058072312558 |
| Transport temperature | Room temperature |
| diameter | 35 mm |
| Color | transparent (clear) |
| Material | Borosilicate glass |
| Volume | 30 ml |
| Form | high shape |
| Height | 50 mm |
| Type closure/ lid | NS lid plug |

