



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



---

**Lenz weighing glass with NS lid stopper, Ø 25 x  
H. 25 mm, ground joint NS 24/12, low form**

**€96.18  
plus VAT &  
Shipping**

**Product Images**

---

## Description

---

### Weighing glasses with NS lid stopper

Low form

With replaceable NS lid stopper, with labeling field. Dimensions (mm) and contents (ml) permanently burned in.

## Additional Information

---

No.	LZ-1522
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3525
EAN	4058072312619
Transport temperature	Room temperature
diameter	25 mm
Color	transparent (clear)
Material	Borosilicate glass
Volume	6 ml
Form	low form
Height	25 mm
Type closure/ lid	NS lid plug

