



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Lenz weighing glass with NS lid stopper, Ø 40 x
H. 25 mm, ground joint NS 40/12, low form**

**€102.42
plus VAT &
Shipping**



Description

Weighing glasses with NS lid stopper

Low form

With replaceable NS lid stopper, with labeling field. Dimensions (mm) and contents (ml) permanently burned in.

Additional Information

No.	LZ-1523
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3526
EAN	4058072312626
Transport temperature	Room temperature
diameter	40 mm
Color	transparent (clear)
Material	Borosilicate glass
Volume	15 ml
Form	low form
Height	25 mm
Type closure/ lid	NS lid plug

