



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Lenz weighing glass with NS lid stopper, Ø 50 x
H. 30 mm, ground joint NS 50/12, low form**

**€105.66
plus VAT &
Shipping**

Product Images

Description

Weighing glasses with NS lid stopper

Low form

With replaceable NS lid stopper, with labeling field. Dimensions (mm) and contents (ml) permanently burned in.

Additional Information

No.	LZ-1525
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3528
EAN	4058072312640
Transport temperature	Room temperature
diameter	50 mm
Color	transparent (clear)
Material	Borosilicate glass
Volume	30 ml
Form	low form
Height	30 mm
Type closure/ lid	NS lid plug

