



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Lenz weighing glass with NS lid stopper, Ø 60 x H.30 mm, ground joint NS 60/12, low form

€181.50
plus VAT &
Shipping

Product Images

Description

Weighing glasses with NS lid stopper

Low form

With replaceable NS lid stopper, with labeling field. Dimensions (mm) and contents (ml) permanently burned in.

Additional Information

No.	LZ-1526
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3529
EAN	4058072312657
Transport temperature	Room temperature
diameter	60 mm
Color	transparent (clear)
Material	Borosilicate glass
Volume	45 ml
Form	low form
Height	30 mm
Type closure/ lid	NS lid plug

