



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**Lenz weighing glass with NS cap, Ø 30 x H. 30 mm, ground joint NS 34/12, low form**

**€141.30  
plus VAT &  
Shipping**

## Product Images



## Description

---

### Weighing glasses with NS cap

With replaceable NS cap, with inscription field.

Dimensions (mm) and contents (ml) permanently burned in.

## Additional Information

---

No.	LZ-1528
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3531
EAN	4058072312671
Transport temperature	Room temperature
diameter	30 mm
Color	transparent (clear)
Material	Borosilicate glass
Volume	12 ml
Form	low form
Height	30 mm
Type closure/ lid	NS cap

