



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Lenz weighing glass with NS cap, Ø 38 x H. 70 mm, ground joint NS 40/12, high form

**€154.32
plus VAT &
Shipping**

Product Images



Description

Weighing glasses with NS cap

With replaceable NS cap, with inscription field.

Dimensions (mm) and contents (ml) permanently burned in.

Additional Information

No.	LZ-1537
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3540
EAN	4058072312763
Transport temperature	Room temperature
diameter	38 mm
Color	transparent (clear)
Material	Borosilicate glass
Volume	75 ml
Form	high shape
Height	70 mm
Type closure/ lid	NS cap

